

FENDERPREP™

INDICATION

FenderPrep™ protects the adjacent tooth and surrounding gingiva on most major preparations such as crowns, veneers and inlays.

FUNCTION

FenderPrep™ is a combined wedge and steel plate. It has a thin "boat-shaped" tip which allows it to be inserted into tight proximal spaces at the same time as it condenses the gingival without damages.

USAGE

FenderPrep™ can be inserted buccally or lingually. The steel plate is perforated thus allowing it to be bent before insertion so that it will lean against the adjacent tooth. This provides protection, better space and vision, and saves time with increased safety. The wedge with its conical profile remains in place during the procedure.

NEW!



Available in:

FenderPrep 36/pkg

6027-49