

Innovative Digital Signal Processing and unique software, developed by MedicNRG, guarantees the most accurate and reliable Apex Locator



ApexNRG[™]XFR

Xtra Fine Resolution Apex Locator

To identify the exact working length of the root canal with a tolerance of 0.1 mm.

Precision at your Fingertips!



Automatic adjustment of the measuring current if the canal condition changes. The ApexNRG-XFR gives a precise and reliable result in any canal condition: **Wet**, **Dry** and **Bleeding** canals.

Xtra Fine Resolution of measurement every 0.25 mm

- User friendly and simple to use
- Automatic calibration
- Automatic filtration of momentary contact with crowns or amalgam fillings
- The unique D.S.P. technology allows a compact hardware design and the minimization of the external dimensions and weight
- Automatic shut down
- Turn on mode by touching the file clip and file or lip hook
- **50 month guarantee**

MEDICNRG[™]
Precise Endo Technology

Distributed by: **JS Dental Manufacturing, Inc.**
PO Box 904, Ridgefield, CT 06877
Tel 203-438-8832 Fax 203-431-8485