

# CALASEPT®

Sterile Calcium Hydroxide Paste... Simplified Delivery System

- Temporary root filling material for conventional root treatment.
- Pulp capping.
- Perforation and fracture repair.
- Stepwise excavation.
- Insulation in deep cavities.

Calcium Hydroxide Root Canal Therapy is the safe and effective way of maximizing the success of root canal procedures. CALASEPT® can be delivered to the root apex in seconds. CALASEPT® is pure calcium hydroxide powder mixed with sterile isotonic saline and packed in self contained syringes. The CALASEPT® System eliminates mixing and clean-up, provides a homogeneous pre-mixed paste consistency and ensures accurate placement, total sterility and full potency with every treatment.

CALASEPT® can be easily removed from the root canal by simple irrigation. CALASEPT® stimulates formation of dentine bridges and hard tissue barrier in roots. CALASEPT® can also be used for stepwise excavation since it has an antibacterial effect on soft carious dentine.

#### 6004u-Large Kit

4 Syringes @ 1.5 ml Calasept Paste  
20 Luer Lock Needles

#### 6002u-Small Kit

2 Syringes @ 1.5 ml Calasept Paste  
10 Luer Lock Needles

#### 6006u-Bulk Kit

20 Syringes @ 1.5 ml Calasept Paste  
100 Luer Lock Needles

#### 6001u-Syringes

5 Syringes @ 1.5 ml Calasept Paste

#### 6003u-Needles

100 Luer Lock Needles

#### 6005u-Cvec Book

#### 6007u-Instructional Video

#### COMPOSITION

100 grams of CALASEPT Paste contains:

- Calcium Hydroxide 41.07g
- Barium Sulphate 8.33g
- Sterile Isotonic Saline Solution 50.60g

pH VALUE 12.4

#### STORAGE

- Store at room temperature.
- Avoid exposure to sunlight.
- Do NOT refrigerate.

#### SHELF LIFE

- 2 years from date of manufacture for unopened cartridge. Approximately 45 days for opened cartridge.

Manufactured by Nordiska Dental AB, Sweden



FL-124 Rev 020107

We Care About How You Practice Dentistry

Brought to you by:  
**JDS DENTAL**  
MANUFACTURING, INC.

Visit us at [www.jsdental.com](http://www.jsdental.com) - P.O. BOX 904 - RIDGEFIELD, CT 06877 - 203-438-8832 - 1-800-284-3368 - FAX 203-431-8485