


SonicTab™ offers you the complete line of cleaning tablets designed to save you money and storage space. SONICTAB TABLETS ELIMINATE HEAVY HANDLING, MIXING AND STORAGE OF LIQUID CONCENTRATES.

SonicTab™

The NEXT Generation of Tablet Cleaning Technology.
The first NON-RESIDUE cleaning tablets.



**One tablet
per gallon
is all it
takes**

**Buy any 3 boxes of
SonicTab, SonicTab TS & SonicTab Evac
and get a FREE Office Cleaning Kit.**
All Purpose Cleaner, Glass Cleaner
and Foaming Hand Soap.
a \$35.85 Value at no Charge
Expires 03/31/08

Professional Strength Cleaning Tablets

For a faster and more effective Ultrasonic Cleaning of your Instruments.



Effective in all water temperatures
Completely water soluble
Functional in wide pH range
No residue to clean up
Aqua color with a hint of mint

- * Anionic compound & solvent break down soils.
- * Non-ionic compound clean metal surfaces without corrosion.
- * Polymer protects and preserves metal.
- * 25 tablets per package (25 gallons).

Non-Corrosive, Non-Toxic, Non-Caustic

For the Cleaning of Dental Appliances, Bridges and Dentures.



Removes: Calculus, Tartar,
Tobacco and Food stains.
Safe for Use on: Dentures,
Bridges Instruments and
Orthodontic appliances.

- * Non-ionic, non-foaming compound and alcohol lifts stains off surfaces.
- * Separate compound whitens surfaces.
- * One tablets makes one pint.
- * 25 tablets per package (25 pints).

For the removal of buffing compounds, food particles & foreign matter.

Keep Your Vacuum System Running without Failures.



Removes and prevents the build up
of debris in the evacuation lines.
Non-ionic, non-foaming compounds
break down soils and keeps
your system clean.

- * Does not cause leaching of mercury from amalgam.
- * Quat prevents the growth of disease causing bacteria from organic and inorganic debris.
- * One tablet makes one application.
- * 25 tablets per package (25 quarts).

Special agents protect and extend the life of your vacuum pump.

Office Cleaning Kit (6 tablets per package with a bottle)



Multipurpose Surface Cleaner
1 tablet makes 24oz bottle.



Glass and Polished Metal Cleaner
1 tablet makes 24oz bottle.



Antibacterial Foaming Hand Soap
1 tablet makes 10oz bottle.